


FICHE TECHNIQUE

| | |
|--------------------|--|
| Référence | 004586 |
| Designation | Plaque 96 puits à fond en verre Corning BioCoat™ - traitement poly - D - Lysine |



Description :

Qty./Pk 1 / Pk
Qty./Cs 10 / Cs
Brand Corning®
Plate Format 96-well half area
Plate Color Clear bottom black
Well Bottom Flat
Well Bottom Color Clear
Well Volume 205 µL
Sterile No
Lids Included Oui
Storage 2°C- 8°C
Surface Coating Poly-D-Lysine