

FICHE TECHNIQUE

Référence	004587
Désignation	Corning® BioCoat® Poly-D-Lysine



Description :

- ✓ High optical quality and scratch-resistant glass
- ✓ Glass bottom thickness of 200 µm is well suited for imaging microscopy
- ✓ Well bottom flatness <50 µm to ensure planarity for imaging devices
- ✓ Low background fluorescence and minimal crosstalk to provide the highest possible optical quality for cell based assays
- ✓ Brand Corning®
- ✓ Plate Format 384-well
- ✓ Plate Color Clear bottom black
- ✓ Well Bottom Flat
- ✓ Well Bottom Color Clear
- ✓ Well Volume 82 µL
- ✓ Sterile No
- ✓ Lids Included Oui
- ✓ Storage 2°C- 8°C
- ✓ Surface Coating Poly-D-Lysine

Unité de vente :

10 unités
