

**Corning Incorporated  
Life Sciences**

**Registered  
ISO 9001:2008**

## Product Description

**Catalog Number:** 6569

**Product Description:** Aluminum Sealing Tape for Thermowell 384 Plate

**Component Materials:**  
Sealing Tape - Aluminum film with an acrylic adhesive.

**Product Dimensions:**  
Length - 5.400 in. Width - 3.250 in.  
Tolerances of Dimensions - +/- .015 in.

**Performance Testing:**  
Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.  
**Visual Attributes:** Visual examination of the product.  
**Packaging:** Inspection for accurate labeling, and correct product configuration.

**Lot Number Designation:**  
8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 2

DOMINIQUE DUTSCHER SAS