

Clearline translucent polystyrene reservoirs

Features

- Designed to accommodate all types of pipettors
- Graduation marks on interior walls
- Made of FDA grade plastic
- Natural polystyrene
- Lot Certified RNase/DNase-Free and Non-Pyrogenic
- Pour-off spouts on each corner
- Sterility validated to 10⁻⁶
- External dimensions : 149 x 69 x 27 mm
- Internal dimensions : 126 x 39 x 27 mm
- 1 reservoir / pack



Unit : 80 Reservoirs