

Corning Incorporated  
Life Sciences

Registered  
ISO 9001:2008

## Product Description

**Catalog Number:** 8508

**Product Description:** Corning® 96 Storage Tube Rack With Locking Lid, Empty, Polypropylene, Not Sterile, 5 Racks/Pack, 2 Packs/Case

**Component Materials:**

Rack - Polypropylene  
Lid - Polypropylene

**Attributes:**

Rack Dimensions - 85.5 x 127.8 mm Base; 24.5 mm Height (Rack Only); 45.0 mm Height (with Lid)  
Hole Pattern - 8 x 12 on 9.0 mm centers

**Heating / Cooling:**

This product is designed for use at temperatures ranging from -100°C (148°F) to 110°C (230°F).

**Sterilization:**

Supplied Not Sterile.

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date; Next 2 digits - Last 2 digits of year of manufacture; Last 3 digits - Batch identification.

Rev No: 1