

Inoculating Loops and Needles

I. APPLICATIONS

1. Microbiology
2. Cell colony isolation and removal

II. FEATURES

1. Materials: polystyrene
2. Gamma radiation sterilized
3. Disposable
4. Evan's Blue Dye Method Certification included

III. LOOP/NEEDLE DESIGN

1. Loops available in 1 μ l and 10 μ l sizes
2. Needles are ideal for the removal of single colonies
3. Ultra-smooth loop surface
4. Needles and Loops are packed in a zip-seal bag

IV. PACKAGING

1. Bulk, 25 pieces per zip-seal bag, 1,000 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
65-0001	1 μ l inoculating loop, blue	1,000 pieces/case	Yes
65-0010	10 μ l inoculating loop, orange	1,000 pieces/case	Yes
65-0002	Inoculating needle, yellow	1,000 pieces/case	Yes

VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
66-0010	Inoculating Loop	200.0	130	66.8	O.D. 3.6	O.D. 6.5	O.D.2.7	O.D. 3.1	O.D.3.8
66-0001	Inoculating Loop	198.3	130	66.8	O.D. 0.8	O.D. 3.0	O.D.2.7	O.D. 3.1	O.D. 3.8
66-0002	Inoculating Needle	198.5	131.7	66.8	O.D. 1.5	O.D. 2.7	O.D.3.1	O.D. 3.8	-



65-0010



65-0001



65-0002