

## Inoculating Loops and Needles

### I. APPLICATIONS

1. Microbiology
2. Cell colony isolation and removal

### II. FEATURES

1. Materials: polystyrene
2. Gamma radiation sterilized
3. Disposable
4. Evan's Blue Dye Method Certification included

### III. LOOP/NEEDLE DESIGN

1. Loops available in 1 $\mu$ l and 10 $\mu$ l sizes
2. Needles are ideal for the removal of single colonies
3. Ultra-smooth loop surface
4. Needles and Loops are packed in a zip-seal bag

### IV. PACKAGING

1. Bulk, 25 pieces per zip-seal bag, 1,000 pieces/case

### V. DESCRIPTION

Item	Description	Packaging	Sterile
65-0001	1 $\mu$ l inoculating loop, blue	25pieces/bag, 40bags/case	Yes
65-0010	10 $\mu$ l inoculating loop, orange	25pieces/bag, 40bags/case	Yes
65-0002	Inoculating needle, yellow	25pieces/bag, 40bags/case	Yes

### VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

### VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.