

Tri-Cornered Beakers

I. APPLICATIONS

1. Transfer of chemicals and other liquids

II. FEATURES

1. Materials: Strong, medical-grade polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. Transparent wall permits easy viewing of tube contents
4. Autoclavable

III. BEAKER DESIGN

1. Unique three-corner design
2. Available in 50ml, 100ml, 250ml, 400ml, 800ml and 1,000ml
3. Molded graduations
4. Clearly visible with either clear or turbid samples

IV. PACKAGING

1. Bulk Packaging, 25 pieces/sleeve and 4 sleeves/case

V. DESCRIPTION

Item	Size	Packaging	Sterile
93-0050	50ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0100	100ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0250	250ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0400	400ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0800	800ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-1000	1,000ml Beaker	25 pieces/sleeve, 4 sleeves/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
93-0050	Beaker	67	58	O.D. 31.1	O.D. 43.1	O.D. 46.9	31	57
93-0100	Beaker	84	72	O.D. 38.7	O.D. 53.7	O.D. 58.4	39	71
93-0250	Beaker	101	88	O.D. 49.1	O.D. 68	O.D. 74	49	88
93-0400	Beaker	114.7	99.9	O.D. 52	O.D. 82	O.D. 84	55.6	110.3
93-0800	Beaker	133.1	115.2	O.D. 65	O.D. 104	O.D.105.1	65	137
93-1000	Beaker	140.6	121.7	O.D. 71.6	O.D.112.7	O.D.114.1	70.9	143

