

Tri-Cornered Beakers

I. APPLICATIONS

1. Transfer of chemicals and other liquids

II. FEATURES

1. Materials: Strong, medical-grade polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. Transparent wall permits easy viewing of tube contents
4. Autoclavable

III. BEAKER DESIGN

1. Unique three-corner design
2. Available in 50ml, 100ml, 250ml, 400ml, 800ml and 1,000ml
3. Molded graduations
4. Clearly visible with either clear or turbid samples

IV. PACKAGING

1. Bulk Packaging, 25 pieces/sleeve and 4 sleeves/case

V. DESCRIPTION

| Item | Size | Packaging | Sterile |
|---------|----------------|----------------------------------|---------|
| 93-0050 | 50ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |
| 93-0100 | 100ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |
| 93-0250 | 250ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |
| 93-0400 | 400ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |
| 93-0800 | 800ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |
| 93-1000 | 1,000ml Beaker | 25 pieces/sleeve, 4 sleeves/case | No |

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

| Cat # | Description | A | B | C | D | E | F | G |
|---------|-------------|-------|-------|-----------|-----------|-----------|------|-------|
| 93-0050 | Beaker | 67 | 58 | O.D. 31.1 | O.D. 43.1 | O.D. 46.9 | 31 | 57 |
| 93-0100 | Beaker | 84 | 72 | O.D. 38.7 | O.D. 53.7 | O.D. 58.4 | 39 | 71 |
| 93-0250 | Beaker | 101 | 88 | O.D. 49.1 | O.D. 68 | O.D. 74 | 49 | 88 |
| 93-0400 | Beaker | 114.7 | 99.9 | O.D. 52 | O.D. 82 | O.D. 84 | 55.6 | 110.3 |
| 93-0800 | Beaker | 133.1 | 115.2 | O.D. 65 | O.D. 104 | O.D.105.1 | 65 | 137 |
| 93-1000 | Beaker | 140.6 | 121.7 | O.D. 71.6 | O.D.112.7 | O.D.114.1 | 70.9 | 143 |

