

Cell Strainers

I. APPLICATIONS

1. Cell culture
2. Filtration of various cells

II. FEATURES

1. Materials: polypropylene frame with nylon mesh
2. Gamma radiation sterilized
3. Non-pyrogenic

III. STRAINER DESIGN

1. Strainers available in 40µm, 70µm, and 100µm nylon mesh
2. Polypropylene frame features a molded tab for easy handling
3. Strainers are color coded for easy identification of mesh size
4. Individually wrapped

IV. PACKAGING

1. Individual Wrapped, 100 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
15-1040-1	40µm Purple strainer	Individually wrapped, 100 pieces/case	Yes
15-1070-1	70µm Orange strainer	Individually wrapped, 100 pieces/case	Yes
15-1100-1	100µm Green strainer	Individually wrapped, 100 pieces/case	Yes

VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
15-1040-1	Cell Strainers	O.D. 35.7	O.D. 25.7	O.D. 20.6	O.D. 7.8	14.4	1.5	2.8	O.D. 25.7
15-1070-1	Cell Strainers	O.D. 35.7	O.D. 25.7	O.D. 20.6	O.D. 7.8	14.4	1.5	2.8	O.D. 25.7
15-1100-1	Cell Strainers	O.D. 35.7	O.D. 25.7	O.D. 20.6	O.D. 7.8	14.4	1.5	2.8	O.D. 25.7

