

## Product Details



### 2.0ml Internal Thread Cryov Traditional

Brand: CryoKING

Material: PP

Sterile: Yes

Temperature: -196°C—121°C

Resistance:

Volume: 2.0ml



2.0ml internal thread, multi codes, sterile O-ring free designed cryogenic vials

### Characteristic

- Made of medical grade polypropylene; can be repeatedly frozen and thawed
- 2.0ml cryogenic vials are available with either internal or external screw threads
- Internal thread vials' caps are with o-rings, which will promote more secure sealing
- DNase & RNase free, endotoxin free, foreign DNA free
- With side barcode and digital code by Side laser printing, easy for information storage

Vials without side barcode also available

- Working temperature: stable from -196°C to 121°C
- Suitable for liquid nitrogen freezing (gas phase)

88-6217	Internal Thread	Yes	Self-Standing	Yes	Yes	Clear	25 Sets/Bag, 40 B.	

DOMINIQUE DUTSCHER SAS