

Cryogenic Vial Closure Color Coders

I. APPLICATIONS

1. Color-coding for cryogenic vials with either internal or external screw threads
2. Marking samples dispensed

II. FEATURES

1. Material: Strong, medical-grade polypropylene, and can be repeatedly frozen and thawed
2. Non sterile
3. Autoclavable to 121°C and freezable to -196°C, and can be stored in vapor phase liquid nitrogen

III. DESIGN

1. Universal screw threads, sealing ring cap
2. Fit cryogenic vials with either internal or external screw threads for quick and convenient classification and visual identification
3. 6 colors: white, red, yellow, blue, green and violet

IV. PACKAGING

Item	Description	Packaging	Sterile
88-6000	Cryogenic Vial Closure Color Coders, White	100 /Bag, 30 Bags/Case	No
88-6001	Cryogenic Vial Closure Color Coders, Red	100 /Bag, 30 Bags/Case	No
88-6002	Cryogenic Vial Closure Color Coders, Green	100 /Bag, 30 Bags/Case	No
88-6003	Cryogenic Vial Closure Color Coders, Blue	100 /Bag, 30 Bags/Case	No
88-6004	Cryogenic Vial Closure Color Coders, Yellow	100 /Bag, 30 Bags/Case	No
88-6005	Cryogenic Vial Closure Color Coders, Violet	100 /Bag, 30 Bags/Case	No
88-6006	Cryogenic Vial Closure Color Coders, Assorted Colors	100 /Bag, 30 Bags/Case, 5 Bags for each color	No

V. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

VI. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.