

## Cryogenic Vial Closure Color Coders

### I. APPLICATIONS

1. Color-coding for cryogenic vials with either internal or external screw threads
2. Marking samples dispensed

### II. FEATURES

1. Material: Strong, medical-grade polypropylene, and can be repeatedly frozen and thawed
2. Non sterile
3. Autoclavable to 121°C and freezable to -196°C, and can be stored in vapor phase liquid nitrogen

### III. DESIGN

1. Universal screw threads, sealing ring cap
2. Fit cryogenic vials with either internal or external screw threads for quick and convenient classification and visual identification
3. 6 colors: white, red, yellow, blue, green and violet

### IV. PACKAGING

Item	Description	Packaging	Sterile
88-6000	Cryogenic Vial Closure Color Coders, White	100 /Bag, 30 Bags/Case	No
88-6001	Cryogenic Vial Closure Color Coders, Red	100 /Bag, 30 Bags/Case	No
88-6002	Cryogenic Vial Closure Color Coders, Green	100 /Bag, 30 Bags/Case	No
88-6003	Cryogenic Vial Closure Color Coders, Blue	100 /Bag, 30 Bags/Case	No
88-6004	Cryogenic Vial Closure Color Coders, Yellow	100 /Bag, 30 Bags/Case	No
88-6005	Cryogenic Vial Closure Color Coders, Violet	100 /Bag, 30 Bags/Case	No
88-6006	Cryogenic Vial Closure Color Coders, Assorted Colors	100 /Bag, 30 Bags/Case, 5 Bags for each color	No

### V. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

### VI. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.