


200 µl (1-330 µl) OneTouch Barrier - Catalog #10350



Catalog #	Description
10350	Natural, Polypropylene, Sterile
Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack
Attributes	



Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free

Material: Polypropylene

Volume Range: 1-330 µl

Grad Lines: 10, 50, 100 µl