


1000 µl (50-1250 µl) OneTouch - Catalog #10360



Catalog #	Description
-----------	-------------

10360	Natural, Polypropylene
-------	------------------------

Packaging Type	Packaging Details
----------------	-------------------

Racked Tips	Non-Sterile, 96 Tips/Rack, 8 Racks/Pack, 5 Packs/Case, Hinged Rack
-------------	--

Attributes
------------



Certified Human DNA Free  
 Certified PCR Inhibitor Free  
 Certified Pyrogen Free

**Volume Range:** 50-1250 µl

**Grad Lines:** 200, 500, 1000 µl