

Digital Hand-held "Pocket" Glucose Refractometer PAL-15S



PAL-15S is suitable to measure the concentration of Glucose water solution.

Specifications

Model	PAL-15S
Cat.No.	4415
Scale	Glucose (g/100g)
Range	0.0 to 50.0% (g/100g)
Accuracy	±0.2%
Temperature compensation	10 to 40° C

Accessories

- PAL-CASE : RE-39409
- JORDAN for PAL(Silicon Cover) : RE-39414
- STRAP : RE-39410
- Battery Case Cover With O-Ring : RE-38163
- USB Battery Charger : RE-34801
- 10% Sucrose (±0.03%) : RE-110010

Patent for Invention Registration No. ZL200310103015.2 (China)

Patent for Design Registration No. ZL03303431.1 (China), 089244 (Taiwan)

Patent Granted in countries around the world.