

Digital Hand-held PAL-19S



PAL-19S is suitable to measure the concentration of Lactose water solution.

Specifications

Model	PAL-19S
Cat.No.	4419
Scale	Lactose (g/100g)
Range	0.0 to 20.0% (g/100g)
Accuracy	±0.2%
Temperature compensation	10 to 40° C

Accessories

- PAL-CASE : RE-39409
- JORDAN for PAL(Silicon Cover) : RE-39414
- STRAP : RE-39410
- Battery Case Cover With O-Ring : RE-38163
- USB Battery Charger : RE-34801
- 10% Sucrose (±0.03%) : RE-110010