

Product Specification Sheet

Item Number: 12180

Description: 2.0 mL Microcentrifuge Tube, Clear, Low Binding

Quantity: 500 Tubes/Bag; 2 Bags/Case

Material: Polypropylene #5

Features:

- RCF Rating of 24,000 x g.
- Graduation lines every .25 mL.
- Frosted area for easy writing on tube.
- Easy close and open cap.
- Available in many colors including blue, green, red, purple, yellow, and amber.

Certifications:

- Non-Pyrogenic, Human DNA Free, Rnase/Dnase Free, BSE/TSE Compliant to EMA/410/01.
- Produced in ISO 13485:2016 Certified Facility.
- Produced in United States of America, compliant to all local and federal laws including EPA and OSHA.

