


## 200 µl (1 - 200 µl) MultiFlex Round - Catalog #13810



Catalog #	Description
13810	Natural, Polypropylene, Sterile
Packaging Type	Packaging Details
Gel Loading Tips	Sterile, Racked, 200 Tips/Pack, 4 Packs/Case, Flex Rack
Attributes	



Certified RNase/DNase Free  
 Certified Human DNA Free  
 Certified PCR Inhibitor Free  
 Certified Pyrogen Free

**Material:** Polypropylene

**Volume Range:** 1-200 µl

**Grad Lines:** 5 µl

**Orifice Outside Diameter:** 0.57 mm