




200 µl (1 - 200 µl) MultiFlex Round - Catalog #13810



Catalog #	Description
-----------	-------------

13810	Natural, Polypropylene, Sterile
-------	---------------------------------

Packaging Type	Packaging Details
----------------	-------------------

Gel Loading Tips	Sterile, Racked, 200 Tips/Pack, 4 Packs/Case, Flex Rack
------------------	---

Attributes



Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free

Material: Polypropylene

Volume Range: 1-200 µl

Grad Lines: 5 µl

Orifice Outside Diameter: 0.57 mm