


20 µl (1-20 µl) Barrier - Catalog #14210



Catalog #	Description
14210	Natural, Polypropylene, Sterile

Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack

**Attributes**



Certified RNase/DNase Free  
 Certified Human DNA Free  
 Certified PCR Inhibitor Free  
 Certified Pyrogen Free

**Material:** Polypropylene  
**Volume Range:** 1-20 µl  
**Grad Lines:** No Graduations

DOMINIQUE DUTSCHER SAS