

TECHNICAL DATASHEET

PART NUMBER

015020CL




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	015020CL
Packaging	ClearLine rack with slip lid
Volume	10 μ l
Graduation	0,5, 10
Length	45,5 mm
Color	Natural
Filter	Yes
Low Retention	No
Sterility	Yes
Autoclavable	No



Sales units :

960 tips