

Thermo Scientific™

Tubes individuels 0.2 ml

Numéro de catalogue: AB0620

Applications liées: PCR



Description

The Thermo Scientific 0.2 mL Individual Tubes with integral “snap shut” caps are suitable for 0.2 mL thermal cycler blocks. They have a 0.25 mL maximum volume when closed. These individual tubes feature thin-wall design and are made from virgin polypropylene. The tubes are cleanroom produced, ensuring that they are free from contaminants that could affect experimental results.

Features of Thermo Scientific 0.2 mL Individual Tubes with Flat Cap, yellow

- Suitable for 0.2 mL thermal cycler blocks
- 0.25 mL maximum tube capacity when capped
- Integral ‘snap shut’ cap

Spécifications

Reaction Speed:	Standard
Type de couvercle.:	Flat
Couleur:	Clear
Volume (métrique):	0.2 mL
À utiliser avec (équipement):	Thermal Cyclers
Type de tube:	Individual