


Thermo Scientific™

Low Profile Tubes and Flat Caps, strips of 8



Catalog number: AB0773

Related applications: PCR

Specifications

Reaction Speed:	Low profile
Type:	Strips of 8; 10 Packs of 12
Cap Type:	Flat
Color:	Clear
Volume (Metric):	0.2 mL
For Use With (Equipment):	Thermal Cyclers
Tube Type:	Strip of 8

Description

Thermo Scientific Strips of 8 Low Profile Tubes with Flat Caps reduces the 'dead space' between the heated lid of the thermal cycler and the sample. This eliminates condensation on the side wall of the tube, preventing reduction in PCR volume and increasing the efficiency of the reaction. Low profile products are especially recommended for use with reaction volumes below 20 μ L.

Features of Thermo Scientific Strips of 8 Low Profile Tubes with Flat Caps

- Suitable for 0.2 mL thermal cycler blocks
- 150 μ L maximum tube capacity when capped
- Low profile to reduce dead space and increase PCR efficiency
- Labelled A–H end tabs (8-strip)
- 10 packs of 12 strips

Applications

- PCR