

Microplate Seal

Low strength adhesive film, transparent, peelable; suitable for short term storage

- ✔ This transparent polyester-based film has a low strength adhesive
- ✔ It is designed as a low-cost sealing option, and useful for temporary storage and as a cover for applications such as centrifugation
- ✔ End tabs allow for easy application and removal
- ✔ This seal is removable without residue on the plate
- ✔ For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)



Seal Integrity
Temperature Range

-20°C
80°C

Peelable
✔

Key Features

- ✔ Non-pierceable
- ✔ Seal integrity range: -20°C to 80°C
- ✔ Peelable
- ✔ Free from DNase, RNase, and human genomic DNA

Use

- ✔ Applications: sample storage (aqueous samples only)
- ✔ Suitable for all plate types

Options

- ✔ Sheet format: 130 x 80 mm

Ordering Information

Code	Details	Quantity
4ti-0510	sheets (130 mm x 80 mm)	100 sheets

DOMINIQUE DUTSCHER