

## Pierce Heat Seal



Seal Integrity  
Temperature Range

-20°C  
120°C

Autoclavable



Pierceable



### Pierce Heat Seal

Pierceable heat sealing foil, high solvent resistance, resealable, suitable for PCR/storage/shipping

Code	Details	Sheets or Rolls/Case
4ti-0530	Roll (610 m x 78 mm, approx. 5,000 seals) <sup>1</sup>	1
4ti-0530/122	Roll (122 m x 78 mm, approx. 1,000 seals) <sup>1</sup>	1
4ti-0530/REMP	Roll (610 m x 78 mm, approx. 5,000 seals) <sup>3</sup>	1
4ti-0530S	Sample roll (5 m x 78 mm) <sup>1</sup>	1
4ti-0532	Roll (500 m x 115 mm, approx. 6,200 seals) <sup>2</sup>	1
4ti-0532/REMP	Roll (500 m x 115 mm, approx. 6,200 seals) <sup>3</sup>	1
4ti-0532S	Sample roll (5 m x 115 mm) <sup>2</sup>	1
4ti-0531	Sheets (125 mm x 78 mm)	100
4ti-0531/GR	Sheets as above, with printed grid reference	100

### Features

- Compatible with PP and PS plates
- Suitable for PCR
- Good solvent resistance - suitable for low temperature and room temperature compound storage in DMSO and organic solvents
- Available also in sheet format with printed grid reference

DOMINIQUE BUREAU BIOTECHER SAS