



Ultimate

high-performance PCR consumables

Benefits

- Polycarbonate frame for warp-resistant thermal cycling
- Enhanced mechanical stability for robotic handling
- Thin-walled wells for optimal heat transfer
- Well shape optimized for maximum sample recovery
- Multiple frame color options, all available with both clear (for PCR) and white colored wells (for qPCR)

Thermo Scientific™ Armadillo™ PCR plates combine the rigidity of a polycarbonate frame with thin-walled polypropylene wells to provide superior thermal cycling performance under all conditions without warping. Armadillo plates are available in 96- and 384-well formats and multiple colors, and they can be ordered with a standard 128 barcode or custom barcoding. All formats are available with clear or white wells, with white giving superior qPCR results. The specially designed warp-resistant frame and multiple format options make Armadillo PCR plates the ultimate choice for high-throughput PCR and automated handling.


**Now
available**
in low &
standard profile
semi-skirted
formats



Thermo
SCIENTIFIC

Instrument Compatibility

Armadillo PCR Plates	Instruments																						
	PCR										qPCR												
	Thermo Scientific [™] Arktik [™] Thermal Cycler	Thermo Scientific [™] PCR Express/ Px2/ PxE/ MBS	Agilent Gradient Cycler	Agilent SureCycler 8800	Biometra TGradient/ TAdvanced	Biometra TRobot/ TProfessional	Bio-Rad T100 [™]	Bio-Rad C1000 [™] / S1000 [™] / PTC-2(XX)/ Dyad/ Tetrad [™]	Eppendorf Mastercycler [®] ep gradient, Pro and Nexus	Gene Technologies GS1/ GS4/ GSX	MWG Primus 96 / 384	Techné [®] Flexigene/ TC-412/ TC-4000	Techné [®] Genius Touchgene/ TC-512/ TC-5000	Techné [®] TC-PLUS/ Prime/ Prime Elite	Bio-Rad CFX96 [™] / CFX384 [™]	Bio-Rad iCycler/ IQ [™] 4/ IQ [™] 5/ MyIQ [™] / MyIQ2 [™]	Bio-Rad Opticon/ Opticon [®] 2, Chromo4 [™]	Eppendorf Mastercycler [®] ep realplex	MWG, The Q Lifecycler	Roche LightCycler [®] 480	Roche LightCycler [®] 96	Techné [®] Prime0	
96-Well Skirted	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
96-Well Low Profile Semi-Skirted	●	●			●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
96-Well Semi-Skirted	●	●	●		●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
384-Well	●	●	●		●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●



All plates are available with random barcodes, but we can customize to your specifications

Fully skirted Armadillo PCR plates

Part No.	Description	Quantity
AB-2396	Armadillo 96-well, clear wells	25 Plates
AB-3396	Armadillo 96-well, white wells	25 Plates
BC-2396	Armadillo 96-well, clear wells, barcoded	25 Plates
BC-3396	Armadillo 96-well, white wells, barcoded	25 Plates
AB-2384	Armadillo 384-well, clear wells	50 Plates
AB-3384	Armadillo 384-well, white wells	50 Plates
BC-2384	Armadillo 384-well, clear wells, barcoded	50 Plates
BC-3384	Armadillo 384-well, white wells, barcoded	50 Plates

Semi-skirted Armadillo PCR plates

Part No.	Description	Quantity
AB-2496	Armadillo 96-well, low profile semi-skirted, clear wells	25 Plates
AB-3496	Armadillo 96-well, low profile semi-skirted, white wells	25 Plates
BC-2496	Armadillo 96-well, low profile semi-skirted, clear wells, barcoded	25 Plates
BC-3496	Armadillo 96-well, low profile semi-skirted, white wells, barcoded	25 Plates
AB-2596	Armadillo 96-well, semi-skirted, clear wells	25 Plates
AB-3596	Armadillo 96-well, semi-skirted, white wells	25 Plates
BC-2596	Armadillo 96-well, semi-skirted, clear wells, barcoded	25 Plates
BC-3596	Armadillo 96-well, semi-skirted, white wells, barcoded	25 Plates

*To select a color please suffix the Part. No. with either /B for blue, /G for green, /O for orange, /R for red or /Y for yellow.

© 2014 Thermo Fisher Scientific Inc. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Indicated trademarks are the property of their respective owners. This product is sold under license from Eppendorf AG.

• For more information, visit:
thermoscientific.com/armadillo

Europe

Customer Service
 cs.molbio.eu@thermofisher.com

Technical Support
 ts.molbio.eu@thermofisher.com

Tel 00800 222 00 888
 Fax 00800 222 00 889

MB1469

United States

Customer Service
 cs.molbio@thermofisher.com

Technical Support
 ts.molbio@thermofisher.com

Tel 877 661 8841
 Fax 800 854 5395

Canada

Customer Service
 cs.molbio@thermofisher.com

Technical Support
 ts.molbio@thermofisher.com

Tel 800 340 9026
 Fax 800 472 8322

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand