

ClearLine® Petri Dishes



Features

- Petri dishes Ø 90 without vent for microbiological analyzes
- Material: Transparent crystal polystyrene
- Bottom Ø: 87.5 mm
- Lid Ø: 91 mm
- Total height: 14 mm
- Flatness: < 1 mm
- Operating temperature: -40°C / +65°C
- Casting temperature: < 60°C

Packaging

- Number of boxes / carton: 500
- Number of boxes / bag: 20
- Storage condition: Room temperature

Quality

- Aseptic
- Traceability: Barcode label, manufacturing batch
- Compliance with standards: ISO 24998
- Absence of heavy metals - directive CE/94/62
- Material food - regulation (EC) n° 1935/2004
- Single use

Units :

500 Petri dishes