


ClearLine® Petri Dishes

Description :

Polystyrene Petri dishes for microbiology, bacteriology

Features :

- Crystal polystyrene, superior clarity
- Manufacturing under aseptic conditions Class ISO 7
- Sterile, irradiated by E-beam, 10kgy
- Sub-packaged in bags
- CE marking, in accordance with the Directive 98/79/EC
- Absence of heavy metals, in accordance with the directive EC/94/62
- Suitable for food contact, according to EC regulation n°1935/2004
- Casting temperature 60°C



Model	Ø 90mm
Vents	Without
Internal dimensions (Ø x h)	88 x 12.7 mm
External dimensions (Ø x h)	91.3 x 13.8 mm
Poids	13,3 g
Units / bag	20

Units :

500 Petri dishes