


Product description

=====

Item-No.: 1020029
Item: Foam swabs with big foam head, wooden applicator, non-sterile



Materials:

Stick: wooden
Head: cotton & foam PPI 80
Fixation: heat

Description: total length with head: A: 150 +/- 3mm
stick: B: 146 +/- 1mm length
D: 2.2 +/- 1mm diameter
head: E: 10 +/- 1mm width
C: 22 +/- 3mm length
F: 5 – 8mm thickness



Intended Use: wound care application, for single use only, non-sterile, only for use through professional medical staff

Shelf life: 5 years

Packaging: 25 pcs./bag, 40 bags/carton, 144 cartons / pallet

Storage: keep in original closed packaging in normal conditions

Instruction for use: no instruction for use needed as the product can be used safely without such instruction (as per annex I, §13 MDD 93/42/EEC). Use only through professional medical staff.

Classification: medical device class I, non-sterile

Norms: DIN EN 13485, DIN EN 980, DIN EN 14971, DIN EN 10993, EN 14079, MPG (German Medical product law), MDD 93/42/EEC, MDD 2007/47/EC