


Product description

Item-No.: 1020011
Item: Foam-Swabs, flexible head, PS-stick, non-sterile



Materials:

Stick: PS (polystyrene)
Head: Foam PPI 80
Fixation: through glue

Description: Total length with head: **A**: 135 +/- 3 mm
Stick: **B**: 130 +/- 1 mm length
D: 5 +/- 1 mm width
G: 2,5 +/- 0,5 mm
Head: **E**: 13 +/- 1 mm width
C: 21 +/- 3 mm length
F: 5 - 8 mm thickness



Intended Use: wound care application, for single use only, non-sterile, only for use through professional medical staff

Shelf life: 5 years

Packaging: 25 pcs./bag, 40 bags/carton, 144 cartons / pallet

Storage: keep in original closed packaging in normal conditions

Instruction for use: no instruction for use needed as the product can be used safely without such instruction (as per annex I, §13 MDD 93/42/EEC). Use only through professional medical staff.

Classification: medical device class I, non-sterile

Norms: DIN EN 13485, DIN EN 980, DIN EN 14971, DIN EN 10993, EN 14079, MPG (German Medical product law), MDD 93/42/EEC, MDD 2007/47/EC