



**Export and Production**  
**- since 1949 -**

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

## **FICHE TECHNIQUE**

---

Item: Steristoppers® No. 8P  
Item No.: 1010801

Material: 2-ply cellulose made of 100% chlorine-free cellulose  
Sinking time: Max. 5 sec.  
Ability of water holding: Min. 13 g/g

Glue: Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it

Shape: Mould twisting process

Frequency of use: For single use

Measuring tolerances: Due to softness of the material

Sterilization: Autoclavable or dry heat up to 180-200°C

Aspiration: Steristoppers® are air permeable

Experience: Steristoppers® have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction

Patents: German Patent No. 1467937 and foreign patents registered.

HS-code: 4823 7090

Shelf Life: N/A

**\*\* This is a computer generated copy of the original document. Proof on file. No signature required. \*\***

**www.herenz.de**

Phone +49-040-739-20-40  
Fax +49-040-739 20-439  
Email hhh@herenz.de  
ILN 40 49434 00000 8

**Managing Directors**

Dieter Cyll, Antje Herenz-Cyll  
Comm. Rg. Hamburg HRB 48884  
Tax No. 44/780/02888  
V.A.T. No. DE 118 665 646

**Postbank Hamburg**

Sort Code 200 100 20  
Account No. 6846 46 200  
IBAN DE97200100200684646200  
Swift-Code PBNKDEFF

**Deutsche Bank AG**

Sort Code 200 700 00  
Account No. 55 00 251  
IBAN DE82200700000550025100  
Swift-Code DEUTDEHH

**Hamburger Sparkasse**

Sort Code 200 505 50  
Account No. 1085 243 002  
IBAN DE87200505501085243002  
Swift-Code HASPDEHH