



**Export and Production
- since 1949 -**

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

Product description

Item: Steristoppers®

Material: 2-ply cellulose made of 100% chlorine-free cellulose
Sinking time: max. 5 sec.
Ability of water holding: min. 13 g/g

Glue: Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it

Shape: Mould twisting process

Frequency of use: for single use

Measuring tolerances: due to softness of the material

Sterilization: autoclavable or dry heat up to 180-200°C

Aspiration: Steristoppers® are air permeable

Experience: Steristoppers® have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction

Patents: German Patent No. 1467937 and foreign patents registered.

Manufacturing Lot No.: printed on each carton: week / year

Shelf Life: N/A

www.herenz.de

Phone +49-040-739-20-40
Fax +49-040-730 41 48
Email hhh@herenz.de
ILN 40 49434 00000 8
Steristopfen-e

Managing Directors

Dieter Cyll, Antje Herenz-Cyll
Commercial Register
Hamburg HRB 48884
Tax No. 44/729/00248
V.A.T. No. DE 118 665 646

Postbank Hamburg

Sort Code 200 100 20
Account No. 6846 46 200
IBAN DE97200100200684646200
Swift-Code PBNKDEFF

Deutsche Bank AG

Sort Code 200 700 00
Account No. 55 00 251
IBAN DE82200700000550025100
Swift-Code DEUTDEHH

Hamburger Sparkasse

Sort Code 200 505 50
Account No. 1085 243 002
IBAN DE87200505501085243002
Swift-Code HASPDEHH