



Export and Production
- since 1949 -

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

Product description

=====

Item: Steristoppers[®], cellulose stoppers

Material: 2-ply cellulose made of 100% elementary-chlorine-free cellulose
Sinking time: max. 5 sec. (DAB 2010)
Ability of water holding: min. 13 g/g (DAB 2010)

Glue: Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it

Shape: Mould twisting process

Frequency of use: for single use

Measuring tolerances: due to softness of the material

Sterilization: autoclavable or dry heat up to 180-200°C

Aspiration: Steristoppers[®] are air permeable

Experience: Steristoppers[®] have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction

Patents: German Patent No. 1467937 and foreign patents registered.

Manufacturing Lot No.: printed on each carton: week / year

Shelf Life: N/A