



**Export and Production**  
**- since 1949 -**

Heinz Herenz Medizinalbedarf GmbH • Rudorffweg 10 • 21031 Hamburg, Germany

## **Product description**

=====

**Item:** Steristoppers<sup>®</sup>, cellulose stoppers

**Material:** 2-ply cellulose made of 100% elementary-chlorine-free cellulose  
**Sinking time:** max. 5 sec. (DAB 2010)  
**Ability of water holding:** min. 13 g/g (DAB 2010)

**Glue:** Glue on basis of methyl-cellulose sprayed on the cellulose while twisting it

**Shape:** Mould twisting process

**Frequency of use:** for single use

**Measuring tolerances:** due to softness of the material

**Sterilization:** autoclavable or dry heat up to 180-200°C

**Aspiration:** Steristoppers<sup>®</sup> are air permeable

**Experience:** Steristoppers<sup>®</sup> have been used to seal tissue cultures for many years in bacteriological institutes all over the world to entire satisfaction

**Patents:** German Patent No. 1467937 and foreign patents registered.

**Manufacturing Lot No.:** printed on each carton: week / year

**Shelf Life:** N/A