



## Dry Swabs without medium

### Dry swabs wooden / PL

Sterile, without medium, plastic tube approx. 150 x 12 mm, **wooden or plastic applicator**, plastic stopper, label.

PU = 100 pcs/box, 3,000 pcs/carton

Item	Description	MOQ
1020001	wooden stick, cotton head	100 pcs
1020002	PS-stick, viscose/rayon head	100 pcs

### Dry swabs Alu

Sterile, without medium, plastic tube approx. 150 x 12 mm, **aluminium applicator**, plastic stopper, label.

PU = 100 pcs/box, 3,000 pcs/carton

Item	Description	MOQ
1020025	aluminium stick, cotton head	100 pcs
1020023	aluminium stick, viscose/rayon head	3,000 pcs

### Dry swabs wooden / PL - peel-pack

Sterile, without medium, plastic tube approx. 150 x 12 mm, **wooden or plastic applicator**, plastic stopper, label, single packed in peel-pack. PU = 150 pcs/box, 900 pcs/carton

Item	Description	MOQ
1020003	wooden stick, cotton head	150 pcs
1020070	PS-stick, viscose/rayon head	900 pcs

### Dry swabs Alu - peel-pack

Sterile, without medium, plastic tube approx. 150 x 12 mm, **aluminium applicator**, plastic stopper, label, single packed in peel-pack. 150 pcs/box, 900 pcs/carton

Item	Description	MOQ
1020071	aluminium stick, cotton head	900 pcs

### Forensic swabs - DNA examined swabs

ETO treated, DNA free, DNASE, RNASE free, **wooden or plastic applicator**. PU = 100 pcs/box, 1,000 pcs/carton

Item	Description	MOQ
1020055	wooden stick, cotton head	100 pcs
1020072	plastic stick, viscose/rayon head	10,000 pcs