


Datasheet

Cylindrical brush

684510-3



Technical details

Product designation	cyl. twisted brush 300/110/30, bent over, with eyelet 13 mm, simple wire handle, galvanized wire 1.6 mm, Fibre
Article no.	684510-3
Total length	300 mm
Brush filling length	110 mm
Brush diameter	30 mm
Brush material	Fibre
Eyelet	13 mm
Version	top bent over
Wire diameter	1.6 mm
Wire material	galvanized wire according to EN ISO 16120-2
Net weight	14.4 g
Packaging unit	-
Country of origin	Germany
Customs tariff number	96039091