



Distributeur mural pour gants 1 boîte
 Dispenser for gloves - wall mounted - 1 box.
 Dispensor mural para guantes - 1 caja

Informations produits / Product information / Informaciones sobre el producto



Référence/ Reference/ Referencia	957210
Matériau / Material / Material	Inox 304 / Stainless steel 304 / Acero inoxidable 304
Code TARIC/ TARIC code / Código TARIC	82159910
Origine/ Origin / Origen	France
EAN 13	3435629572107
Poids / Weight / Peso	0,9 kg
Largeur / Width / Largo	135
Hauteur /Height / Altura	270
Profondeur / Depth / Profundidad	93
Conditionnement / Packging / Condicionado	1
Consommables compatibles / Disposable suitable / Desechables compatibles	Gants de protection jetable / Protective disposable gloves / Guantes de proteccion desechables Nos ref / Our ref / Nustras ref : 957602 - 957603 - 957604 - 657605
Capacité du produit / Product capacity / Capacidad del producto	1 boîte de 100 gants / 1 boxes of 100 gloves / 1 caja de 100 guantes
Fournis avec / Supply with / Suministrado con	Vis et cheville / plugs and screw / tacos y tornillas

Informations colis / Parcel information / Informaciones sobre el paquetes

Poids / Weight / Peso	0,9
Largeur / Width / Largo	145
Hauteur /Height / Altura	275
Profondeur / Depth / Profundidad	100
Volume / Volume / Volumen	0,004 m ³