



Distributeur inox mural pour gants PE
 Wall hanging dispenser for PE gloves in stainless steel.
 Dispensador mural para guantes PE in acero inoxidable

Informations produits / Product information / Informaciones sobre el producto



Référence/ Reference/ Referencia	957600
Matériau / Material / Material	Inox 304
Code TARIC/ TARIC code / Código TARIC	82159910
Origine/ Origin / Origen	France
EAN 13	3435629576006
Poids / Weight / Peso	0,63kg
Largeur / Width / Largo	200
Hauteur/Height / Altura	170
Profondeur / Depth / Profundidad	50
Conditionnement / Packging / Condicionado	2
Consommables compatibles / Disposable suitable / Desechables compatibles	Gants PE sachets de 100 / Gloves PE pack of 100 / Guantes PE bols de 100 Notre ref / Our ref / Nuestra ref : 957601
Capacité du produit / Product capacity / Capacidad del producto	3 - 4 sachets de 100 gants PE / 3- 4 packs of 100 gloves/ 3 - 4 bols de 100 guantes
Fournis avec / Supply with / Suministrado con	Pictogramme / Pictogramm / Pictograma Vis et cheville / plugs and screw / tacos y tornillas Clé / Key / Llave

Informations colis / Parcel information / Informaciones sobre el paquetes

Poids / Weight / Peso	1,69
Largeur / Width / Largo	210
Hauteur/Height / Altura	180
Profondeur / Depth / Profundidad	125
Volume / Volume / Volumen	0,005 m3