





BOX taille L pour casque antibruit, masque respiratoire
 BOX L size for ear defenders, respiratory protection mask
 BOX tamaño L para orejas, mascararas respiratoria

Informations produits / Product information / Informaciones sobre el producto



Référence/ Reference/ Referencia	957952
Matériau / Material / Material	Acier electrozingué + peinture polyester / Galvanized steel + polyester paint / Acero galvanizado + pintura polyester
Code TARIC/ TARIC code /Codigo TARIC	82159910
Origine/ Origin / Origen	France
EAN 13	3435629579526
Poids / Weight / Peso	2,21kg
Largeur / Width / Largo	197
Hauteur /Height / Altura	230
Profondeur / Depth / Profundidad	138
Conditionnement / Packaging / Condicionado	1
Consommables compatibles / Compatible products / Desechables compatibles	Casque anti-bruit, masque respiratoire/ Ear defenders, respiratory protection mask/ Orejas, mascararas respiratorias
Capacité du produit / Product capacity / Capacidad del producto	1 casque anti-bruit ou masque respiratoire/ 1 ear defenders or respiratory mask 1 orejas o mascararas respiratorias
Fournis avec / Supply with / Suministrado con	Pictogramme / Pictogram / Pictograma   Vis et cheville / plugs and screw / tacos y tornillas

Informations colis / Parcel information / Informaciones sobre el paquetes

Poids / Weight / Peso	2,22
Largeur / Width / Largo	210
Hauteur /Height / Altura	260
Profondeur / Depth /Profundidad	175
Volume / Volume / Volumen	0,009m3

Jurine SARL
 95 Rue de La Roche du Geai
 42000 Saint-Etienne
 e-mail : info@jurine.com