




10 µl (0.1 - 10 µl) Micro G - Catalog #23570



Catalog #	Description
23570	Natural, Polypropylene, Sterile
Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 1 Packs/Case, Cell Rack
Attributes	



Volume Range: 0.1-10 µl  
 Grad Lines: 2, 10 µl

Certified Human DNA Free  
 Certified PCR Inhibitor Free  
 Certified Pyrogen Free

DOMINIQUE DUTSCHER