




200 µl (1-200 µl) Standardization - Catalog #26430T



Catalog #	Description
26430T	Natural, Polypropylene, Low Binding, Sterile

Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, HNG Rack

Attributes



Volume Range: 1-200 µl
Grad Lines: 10, 50, 100 µl

Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free