

TECHNICAL DATASHEET

PART NUMBER

027120CL




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	027120CL
Packaging	ClearLine rack with slip lid
Volume	1000 μ l
Graduation	No
Length	72,9 mm
Color	Blue
Filter	No
Low Retention	No
Sterility	Yes
Autoclavable	Yes



Sales units :

1000 tips