

TECHNICAL DATA SHEET

Part Number	027210
Description	Microtube 1.5 ml



Benefits :

- ✓ Material: Polypropylene
- ✓ Colour: natural
- ✓ Bottom: round
- ✓ Skirt: no
- ✓ Marking area: yes
- ✓ Cap: snap cap, designed for repeated opening and closing
- ✓ Low protein adsorption
- ✓ Non sterile
- ✓ Free of DNases, RNases, endotoxins, human DNA and PCR inhibitors
- ✓ Autoclavable

Technical specifications :

- ✓ Graduations: yes
- ✓ Centrifugeable: up to 24,000 g
- ✓ Volume: 1.5ml
- ✓ Dimensions (inner Ø x h without cap): 9.19 x 38.98 mm

Sales unit :

250 Microtubes
