

TECHNICAL DATA SHEET

Part Nbr	028064
Designation	Nalgene plastic Dewar flask



Description :

- ✓ Store samples in liquid nitrogen, dry ice-solvent mixtures or ice water for the short term with Dewar flasks
- ✓ High density polyethylene
- ✓ With high density polyethylene lid
- ✓ Polyethylene coated handle
- ✓ Double wall filled with CFC-free polyurethane foam
- ✓ Use T° : -196°C to +100°C
- ✓ Volume : 2 L
- ✓ Dimensions :
 - Int. Ø 121 mm
 - Working high 225 mm
 - Total high 260 mm

Sales unit :

1 Nalgene plastic Dewar flask