

## TECHNICAL DATASHEET

## PART NUMBER

**028200CL**




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



## ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	028200CL
Packaging	ClearLine rack with slip lid
Volume	10 $\mu$ l
Graduation	2, 10
Length	28,8 mm
Color	Natural
Filter	Yes
Low Retention	No
Sterility	Yes
Autoclavable	No



### Sales units :

960 tips