


200 µl (1 - 200 µl) MultiFlex Round - Catalog #28480



Catalog #	Description
-----------	-------------

28480	Natural, Polypropylene
-------	------------------------

Packaging Type	Packaging Details
----------------	-------------------

Gel Loading Tips	Non-Sterile, Bulk Bag, 1,000 Tips/Bag, 10 Bags/Case , Bulk Pack
------------------	---

Attributes



Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free

Material: Polypropylene

Volume Range: 1-200 µl

Grad Lines: 5 µl

Orifice Outside Diameter: 0.57 mm

DOMINIQUE DUTSCHER SAS