


CAT. NO.	124-0020/0045	
RECEIVER CAP. (mL) ...	115	
PORE SIZE (µm)	0.2/0.45	
MEMBRANE DIA (mm) ...	50	
	INCH	MM
A (UPPER I. D.)	3.20	81.3
B (HEIGHT)	3.77	95.76
C (BASE DIA.)	2.87	72.9
D (OVERALL)	5.22	132.6
NOM. WT. (gms)	62	
MATERIALS		
UPPER/SUPPORT PLATE	PS	
COVER	PS	
RECEIVER	PS	
ADAPTOR	WHITE LOPE	
VENT PLUG	CA	
MEMBRANE	PES PLAIN WHITE	

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponoreno Creek Drive
Rochester, NY 14625-2385, USA

NAME	PB	APP'D	JMC
DATE	JMC	11/29/07	

REV.	1107
REV.	1107

PSS 00-0124-00XX

115mL STERILIZATION
PES FILTER UNITS