


CAT. NO.	125-0020/0045/0080	126-0020/0045/0080	127-0020/0045/0080
PORE SIZE (µm)	0.2/0.45/0.8	0.2/0.45/0.8	0.2/0.45/0.8
UPPER CAP. (mL)	150	250	500
RECEIVER CAP. (mL) ..	150	250	1000
MEMBRANE COLOR	WHITE	WHITE	WHITE
MEMBRANE GRID COLOR (-0020/-0045/-0080) ..	NONE/GREEN/BLACK	NONE/GREEN/BLACK	NONE/GREEN/BLACK
MEMBRANE DIA (mm) ..	50	50	75
	INCHES	MM	INCHES
A (UPPER I.D.)	3.15	80.0	3.80
B HEIGHT	6.15	156.2	11.4
C DIA (BASE)	3.50	88.9	4.48
D OVERALL	3.7	94	4.2
NDM. WT. (gms)	78	95	205

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	GREEN PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LDPE
GASKET	FOAMED LDPE
MEMBRANE	CN
VENT PLUG	CA
PRE-FILTER	GLASS FIBER (127-XXXX ONLY)

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancromo Creek Drive
Rochester, NY 14625-2385, USA

DATE	REV
06/10	0699
DESCRIPTION:	
HC	
PVC	
APP'D	
6/25/99	
PSS 00-12X-XXXX	

STERILIZATION FILTER UNITS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.