


CAT. NO.	155-0020/-0045	157-0020/-0045	156-4020/-4045	158-0020/-0045
PORE SIZE (µm)	0.2/0.45	0.2/0.45	0.2/0.45	0.2/0.45
UPPER CAP. (nL)	150	250	500	500
RECEIVER CAP. (nL)	150	250	500	1000
MEMBRANE COLOR	WHITE	WHITE	WHITE	WHITE
MEMBRANE DIA (mm) ..	50	50	75	75
	INCHES	MM	INCHES	MM
A (UPPER I.D.)	3.15	80.0	3.80	96.5
B HEIGHT	6.15	156.2	9.8	249
C DIA (BASE)	3.50	88.9	3.86	98.0
D OVERALL	3.7	94	4.0	102
NDM. WT. (gms)	78	95	155	205

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	YELLOW PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	SFCA
VENT PLUG	CA
PRE-FILTER	GLASS FIBER (156-XXXX & 158-XXXX ONLY)

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancromo Creek Drive
Rochester, NY 14625-2385, USA

NAME	JG/MG	APP'D	PVC
DATE	PVC	6/25/99	
DESCRIPTION:	PSS 00-15X-XXXX-1		

REV 0699

TISSUE CULTURE FILTER UNITS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.